

Compression Levels

Clinical Indications

Circular Knit
10-15 mmHg
15-20 mmHg

- aching/fatigued legs
- prophylaxis during pregnancy
- mild edema of ankle and foot
- mild varicosities
- post sclerotherapy

Circular Knit
20-30 mmHg
Flat Knit
18-21 mmHg

- aching/fatigued legs
- prophylaxis during pregnancy
- hereditary tendency toward varicose veins
- mild lymphedema
- hypertrophic scar treatment
- prevention/recurrence of venous ulcers
- mild venous insufficiency
- moderate varicosities
- post sclerotherapy
- burn scar management
- DVT prevention

Circular Knit
30-40 mmHg
Flat Knit
23-32 mmHg

- chronic venous insufficiency
- severe varicosities
- post surgical
- post-phlebectomy
- burn scar management
- moderate lymphedema
- prevention/recurrence of venous ulcers
- management of active venous ulcers
- moderate and post traumatic edema
- post sclerotherapy
- pronounced varicosities during pregnancy
- orthostatic hypotension
- prevention of DVT/post thrombotic syndrome

Circular Knit
40-50 mmHg
Flat Knit
34-46 mmHg

- management of active venous ulcers
- severe chronic venous insufficiency
- severe tendency toward edema
- severe lymphedema

Flat Knit
50+ mmHg

- severe post-thrombotic conditions
- severe lymphedema
- elephantiasis

Contra indications

- open or untreated ulcers
- acute thrombophlebitis
- arterial disease
- uncontrolled congestive heart failure
- weeping dermatosis
- phlebothrombosis
- acute dermatitis
- intermittent claudication










For more information or to locate a Juzo dealer, contact us at:

Phone: 1 (800) 222-4999

Fax: 1 (800) 645-2519

www.juzousa.com

support@juzousa.com

Product	Extremity	Stages	Model Number
Juzo Soft 2000 <ul style="list-style-type: none"> extremely soft inconspicuous when worn 		Latency Stage 1	2001 (20-30 mmHg) 2002 (30-40 mmHg)
Juzo Dynamic (Varin) <ul style="list-style-type: none"> soft inside for greater ease in application ideal for sensitive skin durable, opaque knit 		Stage 1 stage 2	3511 (20-30 mmHg) 3512 (30-40 mmHg) 3513 (40-50 mmHg)
Juzo Silver Dynamic <ul style="list-style-type: none"> antimicrobial, anti-static thermodynamic ideal for patients with skin breakdown 		Stage 1 Stage 2	3561 (20-30 mmHg) 3562 (30-40mmHg) 3563 (40-50 mmHg)
Juzo Dynamic (Varin) <ul style="list-style-type: none"> durable, opaque knit wide assortment of styles and sizes 		Stage 1 Stage 2	3511 (20-30 mmHg) 3512 (30-40 mmHg) 3513 (40-50 mmHg)
Juzo Dynamic Full Knit <ul style="list-style-type: none"> durable, opaque knit continued compression in thighs and abdomen 		Stage 2	3512L (30-40 mmHg) 3513L (40-50 mmHg)
Juzo Expert (Helastic) <ul style="list-style-type: none"> dense knit for optimal therapeutic effect can be made in any shape or size 		Stage 2 Stage 3	3021 (18-21 mmHg) 3022 (23-32 mmHg) 3023 (34-46 mmHg) 3024 (50+ mmHg)
Juzo Strong <ul style="list-style-type: none"> micro-massaging weave maximum containment can be made in any shape or size 		Stage 2 Stage 3	3051 (18-21 mmHg) 3052 (23-32 mmHg) 3053 (34-46 mmHg) 3054 (50+mmHg)
Juzo Expert Silver <ul style="list-style-type: none"> Flat knit made with silver Anti-odor, antimicrobial, anti static Machine wash and dry 		Stage 2 Stage 3	3021SV (18-21 mmHg) 3022SV (23-32 mmHg) 3023SV (34-46 mmHg) 3024SV (50+mmHg)
Juzo Silver Strong <ul style="list-style-type: none"> micro-massaging weave maximum containment Eliminates 99% of Staph, Strep and Candida Albicans 		Stage 2 Stage 3	3051SV (18-21 mmHg) 3052SV (23-32 mmHg) 3053SV (34-46 mmHg) 3054SV (50+mmHg)